

RODENSTOCK Instruments



Making analysis
fashionable

Trial Lens Sets

Refraction Accessories

R RODENSTOCK

The Trial Lens Set 266

includes 266 pieces of trial lenses with plastic rims of 38 mm diameter. These lenses perfectly fit for the RODENSTOCK Instruments Trial Frame 10. To complete the Trial Lens Set 266, it comes as a carrying case equipped with a wooden tray and also includes two cross cylinder and a microfiber cloth.

SPECIFICATIONS

Spherical lenses	+/- 20.00D
Cylindrical lenses	+/- 6.00D
Prism lenses	0.50–10.00 PD
Cross cylinder	+/-0.25 // +/-0.50
Accessory lenses	10 pieces
Plastic rim diameter	38 mm
Lens diameter	25 mm




The Trial Frame 10

offers space for a total of 10 trial lenses. You can combine the Trial Frame 10 with the RODENSTOCK Instruments Trial Lens Set 266: the lenses perfectly fit for this trial frame.

SPECIFICATIONS

Axial scale steps	5°
Lens holder load capacity	5 standard trial lenses each on right/left side
Range of PD scale	50–80 mm
Range of nose bridge adjustment	Height: 0–23 mm Angle: 0°–360°
Weight	61 g



The Prism Lens Set 84

is your perfect partner for heterophoria and strabismus measurement. It comes with a variety of 84 prism lenses – four lenses of each prism diopter step. The carrying case of the Prism Lens Set 84 includes a wooden tray and a microfiber cloth.

SPECIFICATIONS

Prism lenses	0.25 – 15.00 PD
Plastic rim diameter	38 mm
Lens diameter	25 mm



The Perimeter Lens Set 68

offers support in determining the field of vision even for people with ametropia. Thanks to the thin metal frame of the lenses, the size of the visual field during perimeter examination is still remained. The 68 lenses are compatible with the lens holders of all perimeter of RODENSTOCK Instruments. A microfiber cloth is included.

SPECIFICATIONS

Spherical lenses	+/- 16.00D
Cylindrical lenses	0.25–6.00D
Metal rim diameter	38 mm
Lens diameter	36 mm

